



AEROCREPE Medium Cotton Crepe Bandage 7.5cm x 4M Wrap/12

\$58.09

In stock
SKU ACM75

 Add to Wish List

 Add to Compare

Medium Crepe Bandage 7.5cm

Details

Size: 7.5cm x 4M **Pricing Unit:** wrap/12 **Carton QTY:** 30 The AEROCREPE Bandage is a high quality cotton crepe bandage, recommended for low support requirements and can be used for the retention of dressings.

Medium Crepe

AEROCREPE medium cotton crepe bandages are of a standard conforming nature which is durable and suitable for dressing retention and medium support requirements.

Easy to Open Wrapper

AEROCREPE medium cotton crepe bandages have a easy open tab built into the wrapper of the bandage. 



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes