



AERODIAGNOSTIC Personal Infrared Ear and Forehead Thermometer

\$199.91

In stock
SKU ADT05

 [Add to Wish List](#)

 [Add to Compare](#)

Ear and Forehead Infrared Thermometer



Details

Pricing Unit: unit/1 **Carton QTY:** 40 The AERODIAGNOSTIC range offers a selection of commonly used diagnostic supplies in the healthcare industry.

Infrared Thermometers

AERODiagnostic Clinical Forehead Infrared Thermometers ensure accurate and stable body temperature readings, without the need for physical contact from a distance of around 5cm.

Simply point the thermometer towards the forehead to read the temperature shown on the LCD screen. The AERODiagnostic Non-Contact Forehead IR Thermometer is specially designed to take the body temperature of a person regardless of room temperature.



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes