



AERODIAGNOSTIC Pistol Grip Infrared Forehead Thermometer

\$167.07

In stock
SKU ADT20

 Add to Wish List

 Add to Compare

AERODIAGNOSTIC Pistol Grip Infrared Forehead Thermometer

Details

Colour: Blue **Pricing Unit:** unit/1 **Carton QTY:** 40

ARTG No.334356

The AERODIAGNOSTIC range offers a selection of commonly used diagnostic supplies in the healthcare industry.

Infrared Thermometers

AERODiagnostic Clinical Forehead Infrared Thermometers ensure accurate and stable body temperature readings, without the need for physical contact from a distance of around 5cm.

Simply point the thermometer towards the forehead to read the temperature shown on the LCD screen.

The AERODiagnostic Non-Contact Forehead IR Thermometer is specially designed to take the body temperature of a person regardless of room temperature.

Specifications:

Product name: Non-contact Thermometer Battery: DC 3V (DC 1.5V AAA 2PCS) Size: 148mm (L)X45mm(W)X90mm(H)
Weight: About 85g (with the battery) Accuracy: Body mode: 0.2(35.0~42.0) Object mode: 0.5(0~100) Measuring range:
Body mode: 32.0~ 42.9 Object mode:0~100 Storage/transportation temperature: -20 to 55 Relative humidity<85%RH
Measuring distant: 3-5cm



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes