



## AEROFORM Conforming Bandage 2.5cm x 4M Wrap/12

**\$16.30**

In stock  
**SKU** AF25

 Add to Wish List

 Add to Compare

Conforming Bandage 2.5cm

### Details

**Size:** 2.5cm x 4M **Pricing Unit:** wrap/12 **Carton QTY:** 60 The AeroForm Conforming bandage is a quality cost-effective bandage with excellent conforming attributes. Our Heavyweight and Premium range offers greater support and compression and are recommended for sprains, strains and snake bites.

### Conforming Specifications

AeroForm Conforming Bandages are of a standard conforming nature which is durable and suitable for dressing retention.

### Easy to Open Wrapper

AeroForm Conforming Bandages have a easy open tab built into the wrapper of the bandage. 



**More Information**

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes