



OPERATOR 5 Series Metal Tough Food & Beverage Kit 42 x 28.5 x 13cm

\$601.97

In stock
SKU AFAK5MF

 Add to Wish List

 Add to Compare


OPERATOR 5 Series Metal Tough Food & Beverage Kit 42 x 28.5 x 13cm

Details

Case: Metal

Size: 42 x 29 x 13cm

Pricing Unit: unit/1

This food-grade kit is designed specifically for the food handling and preparation industries, containing the essential visual detectable plasters and tape. Compliance with workplace first aid regulations is achieved across all states and territories with this broad kit range. Workplace injuries in higher risk areas can be treated effectively with the quality selection of components. **Contents meet Australian Standard AS2675-1983**  This kit is available in a Tough Metal cabinet that can be wall mounted and is key-lockable for security - all our cabinets are available with custom labelling.

This First Aid Kit comes fully stocked and as part of our **Instant Restock™** system you can refill and service your client's kits quickly with refill pouches designed to fit into each case type. They provide quick access to all product through a unique centre slit pocket without the need of unpacking each individual product.

KIT CONTENTS:

| CODE | DESCRIPTION | QUANTITY |
|----------|------------------------------------|----------|
| AV2050-I | AEROPLAST Visual Plaster 75x25mm | 50 pc |
| AB3-I | AEROBURN Burn Gel Sachet 3.5g | 5 pc |
| PTD06-I | AERODIAGNOSTIC Diagnostic Penlight | 1 pc |



| | | |
|------------|---|---------|
| AAS50-I | AEROAID Antiseptic Spray 50ml | 1 pc |
| AFS001-I | AEROSHIELD CPR Face Shield | 1 pc |
| AGNPF02-I | AEROGLOVE Nitrile Examination Gloves | 5 Pairs |
| AGS753S-I | AEROSWAB Gauze Swabs 7.5cm x 7.5cm 3's | 9 Pks |
| AW1000-I | AEROWASH Eye Wash Ampoule 15ml | 8 pc |
| AW8001-I | AEROWIPE Cleansing Wipe | 10 pc |
| ASP001-I | AEROPROBE Splinter Probes 3.7cm | 10 pc |
| APD50S-I | AEROPAD Low Adherent Dressing 5cm x 5cm | 6 pc |
| APD100S-I | AEROPAD Low Adherent Dressing 7.5cm x 10cm | 3 pc |
| APD101S-I | AEROPAD Low Adherent Dressing 10cm x 10cm | 2 pc |
| AF50-I | AEROFORM Conforming Bandage 5cm x 4M | 3 Roll |
| AF75-I | AEROFORM Conforming Bandage 7.5cm x 4M | 3 Roll |
| AFHP100-I | AEROFORM Heavyweight Conforming Bandage 10cm x 4.5m | 2 Roll |
| AV2010-I | AEROPLAST Visual Tape 2.5cm x 5M | 1 Roll |
| AWD13S-I | AEROWOUND BPC Wound Dressing #13 | 2 Roll |
| AWD14S-I | AEROWOUND BPC Wound Dressing #14 | 3 Roll |
| AWD15S-I | AEROWOUND BPC Wound Dressing #15 | 3 Roll |
| ACD1020S-I | AEROWOUND Combine Dressing 10cm x 20cm | 1 pc |
| ABC10-I | AEROBAND Triangular Bandage 110 x 110 x 155cm | 4 pc |
| ATB130-I | AERORESCUE Emergency Rescue Blanket 127 x 180cm | 1 pc |
| AEP1S-I | AEROPAD Eye Pads 5.5cm x 7.7cm | 4 pc |
| AII2500-I | AEROCOOL Instant Ice Pack 80g | 1 pc |
| ND90-I | AEROHAZARD Needle Disposal Unit 90ml Yellow | 1 pc |
| TRYLD1-I | AEROHAZARD Disposable Dressing Tray | 1 pc |
| PPC60-I | AEROHAZARD Portion Cup Plastic 60ml | 1 pc |
| ASC13-I | AEROINSTRUMENTS Scissor 13cm | 1 pc |
| AF13-I | AEROINSTRUMENTS Tweezer 11cm | 1 pc |
| AAB100-I | AEROSUPPLIES Amputated Parts Bags Assorted Sizes | 3 pc |
| ASN096P-I | AEROSUPPLIES Notebook & Pen | 1 pc |
| ASP12-I | AEROPINS Safety Pins | 12 pc |
| AGL001-I | AEROGUIDE First Aid Leaflet | 1 pc |



More Information

| | |
|------------------------|----------|
| MANUFACTURER_MS | Germany |
| MODEL NUMBER_MS | K61SD521 |
| ECOLABEL_MS | yes |