




MODULATOR Purple Dressings & Bandage Module 20 x 10 x 6cm

\$57.13

In stock
SKU AFAKMODD

 Add to Wish List

 Add to Compare

MODULATOR Purple Dressings & Bandage Module 20 x 10 x 6cm



Details

The Modulator Dressing & Bandage Module is a visual and clearly labelled soft pack that can be ordered individually as a replacement module for The Modulator Softpack or The Modulator Metal Cabinet editions.

This module can also be used in addition to any core workplace kit or as a single standing dressing and bandage kit.

Ensures hasslefree servicing and easy re-stocking with individual replacement modules.

Contents meet Australian Standard AS2675-1983



CODE	DESCRIPTION	QUANTITY
AGS753S	AEROSWAB Gauze Swabs 7.5 x 7.5cm 3s	3 pc
APD50S	AEROPAD Low Adherent Dressings 5cm x 5cm	3 pc
APD100S	AEROPAD Adherent Dressings 7.5cm x 10cm	2 pc
APD101S	AEROPAD Adherent Dressings 10cm x 10cm	1 pc



AF50	AEROFORM Conforming Bandage 5cm x 4m	2 Rolls
AWD14S	AEROWOUND BPC Wound Dressing #14	1 pc
AF75	AEROFORM Conforming Bandage 7.5cm x 4m	1 pc



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes