



AEROFILM Waterproof Film Dressing 6 x 7cm Box/20 (GST FREE)

\$19.15

In stock
SKU AFW00620

 [Add to Wish List](#)

 [Add to Compare](#)

AEROFILM Waterproof Film Dressing 6 x 7cm Box/20 (GST FREE)

Details

Size: 6 x 7cm **Pricing Unit:** box/20 **Carton QTY:** 8 AeroFilm is a sterile, waterproof, ultra-thin, transparent film dressing. It is designed to be applied over I.V. devices, clean skin or wounds and provides a breathable anti-bacterial and viral barrier to outside contaminants.

SupaBond adhesive

- Exclusively engineered for unrivalled adhesion, ensures that AeroFilm can be left in place for up to 7 days. This reduces interference with the wound, reducing nursing time and aiding faster recovery.

Waterproof Film

- AeroFilm is impermeable to water and body fluids, protecting the application underneath and allowing the patient to shower and bathe without redressing. It also aids in the prevention of bacterial contamination, and reduces risks of skin damage through interference.

Simple application

- AeroFilm has an integrated frame applicator that provides a simple and effective application system. The frame enables the dressing to be applied crease free and ensures adhesion. AeroFilm transparency enables correct positioning over the wound site.



Adaptable to Awkward Areas

- AeroFilm can be easily applied to the most awkward body contours, due to its expandability, while SupaBond ensures the dressing remains in place.

Contraindications / Precautions

- AEROFILM should not be used on infected wounds unless under the supervision of a healthcare practitioner.
- As with all adhesive products, apply and remove carefully from sensitive or fragile skin.
- Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

This product is not a substitute for direct medical attention from a healthcare practitioner. 



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes