



AEROFORM Size E Adult Legs Elastic Tubular Bandage 8.5cm x 1M

\$8.16

In stock
SKU AGE01

 Add to Wish List

 Add to Compare

AEROFORM Size E Adult Legs Elastic Tubular Bandage
8.75cm x 1M


Details

Size: 8.5cm x 1M **Pricing Unit:** unit/1 **Carton QTY:** 10 AEROFORM Tubular Bandages are an elasticised support bandage suited to treatment of sports injuries, dislocations and sprains. The Size E tubular bandage is suitable for use on adult legs.

Self Supporting

AEROFORM Tubular Bandages are held in place without tying or tapes, resulting in a comfortable long wearing bandage.

Provides Even Support

The AEROFORM Tubular Bandages provides even support across damaged joints and muscles due to its elastic self supporting properties resulting in a comfortable support bandage. Available in 1m lengths or 10m rolls. 



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes