





AEROSWAB Non-Sterile White Gauze Swab 5 x 5cm Pack/100 (8PLY)

\$3.24

In stock
SKU AGS50

-  Add to Wish List
 -  Add to Compare
- 8 ply Gauze Swab



Details

Size: 5 x 5cm **Pricing Unit:** pack/100 **Carton QTY:** 160 The AEROSWAB™ Gauze Swab is highly absorbent cotton gauze swab that when soaked in fluid becomes an ideal item to clean and prepare wounds for dressing.

Cotton Gauze

- AEROSWAB™ cotton gauze has a balanced combination of softness, abrasive texture and absorbency for effectively cleaning wounds.

Multiple Sizes Available

- AEROSWAB™ gauze swabs are available in a number of different sizes to suit any application.





More Information

| | |
|------------------------|----------|
| MANUFACTURER_MS | Germany |
| MODEL NUMBER_MS | K61SD521 |
| ECOLABEL_MS | yes |