



## AEROCASE Red and Black Plastic Tacklebox with 2 Trays 20 x 40 x 23cm

**\$67.67**

In stock  
**SKU** AKP002BR

 [Add to Wish List](#)

 [Add to Compare](#)

Plastic Case Red and Black 2 Tray

### Details

**Size:** 20 x 40 x 23cm **Pricing Unit:** case **Carton QTY:** 6 The AeroCase Range of plastic cases caters to a broad range of first aid kit requirements in the industry.

### Toolbox Style Plastic Cases

These top quality plastic toolboxes are made in Australia and feature lift-out or cantilever trays so that all products are visible and readily accessible.



### More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes