



AEROCASE Red and White Plastic Tacklebox with 2 Trays 20 x 40 x 23cm

\$67.67

In stock
SKU AKP002R

 [Add to Wish List](#)

 [Add to Compare](#)

Plastic Case Red and White 2 Tray

Details

Size: 20 x 40 x 23cm **Pricing Unit:** case **Carton QTY:** 6 The AeroCase Range of plastic cases caters to a broad range of first aid kit requirements in the industry.

Toolbox Style Plastic Cases

These top quality plastic toolboxes are made in Australia and feature lift-out or cantilever trays so that all products are visible and readily accessible.



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes