



AEROSWAB Sterile White Non-Woven Swab 10 x 10cm (Packs of 3)

\$11.43

In stock
SKU ANS103S

 [Add to Wish List](#)

 [Add to Compare](#)

Non-Woven Swab - 3pcs



Details

Size: 10 x 10cm **Pricing Unit:** box/25 **Carton QTY:** 24 The AEROSWAB™ Gauze Swab is highly absorbent cotton gauze swab that when soaked in fluid becomes an ideal item to clean and prepare wounds for dressing.

Cotton Gauze

- AEROSWAB™ cotton gauze has a balanced combination of softness, abrasive texture and absorbency for effectively cleaning wounds.

Multiple Sizes Available

- AEROSWAB™ gauze swabs are available in a number of different sizes to suit any application.

Sterile

- AEROSWAB™ gauze swabs are also available in a sterile form for minimising infection risk when dressing clean wounds





More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes