



AEROPAD Non-Adherent Dressing 10 x 10cm Box/3

\$4.37

In stock
SKU APD10135

 Add to Wish List

 Add to Compare

Low Adherent Dressing Pad Envelope/3



Details

Size: 10 x 10cm **Pricing Unit:** env/3 **Carton QTY:** 96 The AEROPAD Non-Adherent dressing is recommended for treating cuts, grazes and minor burns using specialized SupaSorb technology.

Premium Absorbency

- Premium absorbent inner material for weeping and bleeding wounds

Non-Adherent

- Special non-adherent film on each side of the pad minimizing wound trauma when changing dressings.

Individually packed and sterile

- Dressings are available individually on demand and sterile for wound cleanliness.





More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes