



AEROPAD Non-Adherent Dressing 5 x 5cm Box/50

\$6.40

In stock
SKU APD50S

 Add to Wish List

 Add to Compare

Low Adherent Dressing Pad Box/50



Details

The AEROPAD Non-Adherent Dressing is recommended for treating cuts, grazes and minor burns using specialised SupaSorb technology, sized 5cm x 5cm and supplied in a bulk box of 50 for clinics and workplaces.

Features

- Premium absorbent inner material for weeping and bleeding wounds
- Non-adherent film on each side minimises trauma when changing dressings
- Individually packed and sterile

Specifications

- Size: 5 x 5cm
- Pack: Box of 50



Common uses

Treating minor cuts and grazes, restocking clinics and workplaces with a bulk supply. Need a smaller pack? See our [5x5cm \(Box/5\)](#).

AEROPAD 5x5cm (Box/50) FAQs

What technology does AEROPAD use?

SupaSorb technology, with a premium absorbent inner material for weeping and bleeding wounds.

Is this dressing non-adherent?

Yes - a special non-adherent film on each side minimises wound trauma when changing the dressing.

What size and pack is this?

5 x 5cm, box of 50.

Do you have other sizes?

Yes - see our other AEROPAD Non-Adherent Dressing sizes: [7.5x10cm](#), [10x10cm \(Box/3\)](#), [10x10cm \(Box/50\)](#), [10x20cm](#), [5x5cm \(Box/5\)](#), [5x7.5cm](#)

Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with



any questions that you may have regarding a medical condition or injury.



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes