



AEROPLAST Kinesiology Tape 5cm x 5M Box/12

\$92.95

In stock
SKU APKT50



 Add to Wish List

 Add to Compare

Kinesiology Sports Tape



Details

Size: 5cm x 5M **Pricing Unit:** box/12 **Carton QTY:** 6

Our New Premium Sports Tape Range

- Supabond™Hypoallergenic Adhesive
- High-Quality Performance Based Materials
- Competitive Price
- Latex Free
- Lightweight for Comfort and Performance
- Flexible Kinesiology and Rigid Options
- Trusted on and off the field

AeroPlast Kinesiology Tapes are a highly flexible elastic tape that adheres tightly to skin promoting increased blood flow and relieves pressure on injured muscles. Supabond™Hypoallergenic Adhesive is latex free with a low allergy risk that provides a secure adhesive with excellent conformability for joints in awkward areas.





More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes