



AEROPLAST Rigid Sports Tape 3.8cm x 13.7M Box/12

\$81.07

In stock
SKU APRT38



 Add to Wish List

 Add to Compare

Rigid Sports Tape Medium



Details

Size: 3.8cm x 13.7M **Pricing Unit:** box/12 **Carton QTY:** 6

Our New Premium Sports Tape Range

- Supabond™Hypoallergenic Adhesive
- High-Quality Performance Based Materials
- Competitive Price
- Latex Free
- Lightweight for Comfort and Performance
- Flexible Kinesiology and Rigid Options
- Trusted on and off the field

AeroPlast Kinesiology Tapes are a highly flexible elastic tape that adheres tightly to skin promoting increased blood flow and relieves pressure on injured muscles. Supabond™Hypoallergenic Adhesive is latex free with a low allergy risk that provides a secure adhesive with excellent conformability for joints in awkward areas.

AeroPlast Rigid Sports Tapes are a lightweight and latex free rigid tape, ideal for strapping and supporting strains and sprains. The tape is designed to be hand tearable.





More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes