



## AEROSPLINT Folded Aluminium Foam Splint 90 x 11cm

**\$24.38**

In stock  
**SKU** AR1018F

 Add to Wish List

 Add to Compare

Aluminium Foam Splint Flat



### Details

**Size:** 90 x 11cm **Pricing Unit:** unit/1 **Carton QTY:** 50

### Light Weight Flat Splint

Immobilisation is a critical item of care when responding to a casualty to limit the occurrence of further injury. Quality immobilisation using suitable products is integral to a full recovery of skeletal injuries.

### Aluminium Foam Splints

An aluminium foam splint is a highly adaptive splint that can be rolled folded or formed into almost any application with the use of a quality support bandage. Available in a rolled or flat form, shrink-wrapped for protection.



## More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes