



AEROSPLINT Folded Extra Large Aluminium Foam Splint 91 x 14cm

\$31.37

In stock
SKU AR1091XL

 Add to Wish List

 Add to Compare

AEROSPLINT Folded Extra Large Aluminium Foam Splint 91 x 14cm

Details

AeroSplint Folded Extra Large Aluminium Foam Splint

Quality immobilisation is critical for full recovery of skeletal injuries. The AeroSplint is a highly adaptive splint that can be folded to support a fractured or injured limb. AeroSplint can be folded to better suit the application site, ensuring a secure fit for splinting. Shrink wrapped for protection the Extra Large AeroSplint offers greater support for larger limbs. Other sizes are available.



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes