



HYFIN Vent Chest Seal Pack/2

\$62.98

In stock
SKU ARHYV1

 [Add to Wish List](#)

 [Add to Compare](#)

[HyFin Vent Chest Seal - Twin Pack](#)

Details

The HyFinVent Chest Seal Twin Pack is the industry leader in the treatment of penetrating injuries to the chest. Two separately packaged vented chest seals are offered in one unique package for treatment of both entry/exit or multiple penetrating injuries to the chest. The HyFinVent Chest Seal design provides 3-vented channels that prevent airflow into the chest cavity during inhalation, while allowing air to escape through the vent channels during exhalation. The 3-vent channels allow blood to escape and also provide a backup fail-safe system, as even if two of the three channels become obstructed, the vent will remain fully operational. Advanced adhesive technology provides superior adhesion in the most adverse conditions, including sweaty or hairy casualties. Packaged in rugged, easy-to-open foil pouches, the perforated packaging allows rescuers to open only one dressing at a time as needed. Each chest seal also includes a gauze pad to wipe the wound surface prior to application. Each HyFinVent Chest Seal has a large, pull tab for single-step peel-and-apply application and allows for the burping of the wound if necessary. The clear, transparent backing allows for easy placement over the wound area and conformability to the patients chest.



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes