



## AEROTAPE Transparent Plastic Tape 2.5cm x 5M Box/12

**\$19.87**

In stock  
**SKU** AT25



 Add to Wish List

 Add to Compare

Microperforated Transparent Tape (wrapped rolls)



### Details

**Size:** 2.5cm x 5M **Pricing Unit:** box/12 **Carton QTY:** 24 AEROPLAST Transparent Tape is a flexible and breathable transparent tape suitable for the fixation of medical devices to the skin in hospital and clinical settings.

### FEATURES & BENEFITS

#### Hypoallergenic

AEROPLAST Transparent Tape has a very low irritant risk, and is suitable for use on fragile and delicate skin.

#### Ease of Use

AEROPLAST Transparent Tape is micro-perforated for easy tear-off on any direction and easy adjustment to size.


#### Individually Wrapped

AEROPLAST Transparent Tape comes in boxes of individually wrapped rolls ensuring the tape is kept in optimum condition and ready for use

#### Transparent

AEROPLAST Transparent Tape is transparent allowing for improved aesthetics and the monitoring the cleanliness of the



site without removing the product. 



### More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes