



## AEROWASTE Bio-Waste Vomit Bag 1500ml Bag/4

**\$5.80**

In stock  
SKU AW1514

 Add to Wish List

 Add to Compare

AEROWASTE Bio-Waste Vomit Bag 1500ml (4-Pack) is a compact pack of sealable biohazard waste bags for home first aid kits, vehicle glove-boxes, marine first aid stations and small workplaces. Each bag holds 1500ml with a leak-proof seal and biohazard marking.



### Details

## What is the AEROWASTE Vomit Bag 1500ml (4-Pack)?

The AEROWASTE Bio-Waste Vomit Bag 1500ml 4-Pack is the compact version of our 50-pack bulk box - designed for home first aid kits, vehicle glove-boxes, marine first aid stations and small workplaces. Each bag holds up to 1500ml of vomit or bodily fluids, with a leak-proof seal and biohazard marking. Made in Australia by AERO Healthcare.

## Where the 4-pack fits

- Family home first aid kits.
- Car, 4WD and caravan glove-box first aid kits.
- Marine and boat first aid stations.
- School and childcare backpacks.
- Travel first aid kits for road trips and overseas travel.
- Small workplaces (under 20 staff).
- Sports team first aid kits.



## Key features

- **1500ml capacity per bag** - covers a single typical vomit episode.
- **Leak-proof seal** - twist-tie closure prevents fluid escape during disposal.
- **Biohazard marking** - clearly identifies waste for safe disposal.
- **4 bags per pack** - suits households, vehicles and small workplaces.
- Compact size - fits in glove-boxes, handbags and small first aid cabinets.

## How to use AEROWASTE vomit bags

1. Hand the bag to the patient or hold open in front of them.
2. Allow vomit to collect inside.
3. Twist or gather the top of the bag closed.
4. Tie firmly with the integrated seal.
5. Dispose in general or biohazard waste per local council policy.
6. Wash hands thoroughly after handling.

## Related products

- [AEROWASTE 50-Pack](#) - bulk pack for clinics, ambulances and larger workplaces.
- [AEROWASTE Vomit Bag Dispenser](#) - wall-mount storage.
- For complete kit stocking see [first aid kits](#).

## Frequently asked questions

Why choose the 4-pack instead of the 50-pack?

The AEROWASTE 4-pack is for households, vehicles, boats and small workplaces where a full box of 50 bags would expire before use. The 4-pack fits in a glove-box, console or home first aid cabinet. For ambulances, clinics or workplaces with more than 20 staff, the 50-pack is more economical.

Can I use the 4-pack for car sickness?

Yes. The 1500ml AEROWASTE bags are commonly stocked in family cars and 4WDs for travel sickness, especially with children prone to motion sickness on long drives. The sealable design prevents spills during travel and contains odour until the next service-station bin.



How long do AEROWASTE bags last in storage?

Stored in a clean, dry place at room temperature, AEROWASTE bags have a long shelf life - typically 5 years from manufacture. Check for any damage or seal degradation before use. Replace promptly if punctures or seal weakness are visible.

Are AEROWASTE bags safe for use during pregnancy?

Yes. The bags themselves are inert polyethylene with no chemicals that could harm pregnant women. Many expectant mothers keep AEROWASTE 4-packs in their handbag for morning-sickness episodes outside the home, given the discreet seal-and-dispose design.

Sources: [Safe Work Australia — First aid](#), [healthdirect.gov.au](http://healthdirect.gov.au) — [Vomiting and nausea](#).



### More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes