



4L Biohazard Clinical Waste
Bag 250 x 300mm - Press
Seal, 30um

\$1.43

In stock
SKU BHB2538

 [Add to Wish List](#)

 [Add to Compare](#)

4L Biological Hazard Clinical Waste Bag



Details

Size: 250 x 300mm **Pricing Unit:** bag **Carton QTY:** 1000 Aero offers a wide range of waste disposal containers designed to meet user demands for strength, simplicity and value. This range offers outstanding protection against penetrating sharps and leaking fluids due to its rigid polypropylene construction and wall thickness. AeroHazard Bags are a highly visual and durable waste bag for biohazardous and clinical waste collection, available in a variety of sizes.



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes