




ChoiceMMed Black Finger Pulse Oximeter

\$157.33

In stock
SKU C63

 [Add to Wish List](#)

 [Add to Compare](#)

ChoiceMMed Finger Pulse Oximeter

Details

Pricing Unit: unit/1 **Carton QTY:** 1 The AERODIAGNOSTIC range offers a selection of commonly used diagnostic supplies in the healthcare industry.

Pulse Oximeters

The Pulse oximeter provides exceptional quality and performance in a small, robust design. The devices are commonly used to measure oxygen saturation levels and pulse measurement, to advise if oxygen supply is required and determine heart rate. It is suitable for use on adults as well as paediatrics in a variety of applications like GP surgeries, Clinics, Dental surgeries, Ambulances, Hospitals, Home care through to spot-checking before, during and after exercise.



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes