



CARDI ACT Alarmed AED Cabinet with Strobe Light 42 x 38 x 15.5cm

\$398.98

In stock
SKU CC-50


 [Add to Wish List](#)

 [Add to Compare](#)

CardiAct AED Wall Cabinet with Alarm & Strobe



Details

Size: 42 x 38 x 15cm **Pricing Unit:** unit/1 **Carton QTY:** 1 Suits most conventional AED models, including HeartSine and Lifepak AEDs. 



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes