




CARDIACT Alarmed Outdoor AED Cabinet with Lock 48 x 47 x 31cm

\$542.25

In stock
SKU CC-60

 [Add to Wish List](#)

 [Add to Compare](#)

CARDIACT Alarmed Outdoor AED Cabinet with Lock 48 x 47 x 31cm

Details

AED Alarmed Outdoor Cabinet w/lock 48cm x 47cm x 31cm



Extra Protection

Waterproof sealing restricts moisture from entering the cabinet



Visible

AED signage on sides of cabinet



High Quality

Premium materials with powder coated finish



Secure

Keypad lock restricting unauthorised access



Resistant

Galvanised steel construction protects the cabinet from rust when stored in a wet or salty environment



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes