



Plastic Screening Kit 36 x 25 x 8.5cm

\$259.83

In stock
SKU CPSPK03

 Add to Wish List

 Add to Compare

Plastic Screening Kit 36 x 25 x 8.5cm



Details

To be used at the screening station prior to staff entry. Contains all the essential products needed for setup and clean-up from conducting the daily scanning procedure. This also contains easy to follow instructions as to how screening is to be conducted.



KIT CONTENTS:

ADT10	AERO DIAGNOSTIC Clinical Infrared Forehead Thermometer	1 each
ACG500	AEROCLEANSE Anti-Bacterial Hand Gel 500ml	2 Bottle
AWB7000	AEROWIPE Alcohol Wipes 155 x 140mm	1 tub
AGNPF02	AEROGLOVE Nitrile Glove Powder Free x 2 (Large)	2 5 bag
CM01	AEROMASK Civilian 3Ply Mask W/Ear loops	2 5 each
SG002	AEROSHIELD Safety Spectacles Vortex Clear	2 each
BHB2538	AEROWASTE Clinical Waste Bag 250x300mm - 4L Press 4L Press Seal, 30um	5 bag
CUSTLEAF	Leaflet/Instructions	1 each



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes