



TRAUMASIM Arm Haemorrhage Control Stabilising Bracket

\$287.29

In stock
SKU HCS01A

 Add to Wish List

 Add to Compare

TraumaSim Haemorrhage Control Simulator Arm

Details

This bracket system can be used to secure the TraumaSim Haemorrhage Control Lifecast Arm (HCS01) arm to a bench or table. 2x additional drill holes are available should you wish to permanently secure the bracket and/or arm to a sturdy structure additional bolts & washers are required for this, not supplied.

Includes:

Sturdy & durable plastic bracket

2 x 30mm hex-head bolts

2 x cast iron g-clamps



[Check out the HCS01 TraumaSim Haemorrhage Control Simulator Arm product here](#)



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes