



BRAYDEN Junior CPR Manikin Kit

\$299.24

In stock
SKU IM13-J

 Add to Wish List

 Add to Compare

Brayden CPR Manikin (no illuminators)



Details



Size: Junior

Pricing Unit: unit/1

Carton QTY: 1

Convert a Brayden adult manikin into a junior within seconds using the Brayden Junior Manikin Kit. 

How to Convert an Adult to a Junior:



1
Open up manikin and remove adult spring.

2
Insert spring adaptor into the base of the manikin.

3
Connect junior spring into the spring adaptor.

4
Connect the second spring adaptor to the rib frame.

5
Close the rib section, ensuring prongs are placed into the slots.

6
Replace adult face and body skin with junior face and body skin.



Compatible with:

[Brayden](#), [Brayden Pulse](#) & [Brayden Pro](#)



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes