



BRAYDEN Baby CPR Manikin

\$2,090.93

In stock
Only **290** left
SKU IM17

 Add to Wish List

 Add to Compare

Brayden CPR Manikin (no illuminators)




Details

Size: Baby

Pricing Unit: unit/1

Carton QTY: 1

 The Brayden Baby is the new arrival to our family' range of Brayden manikins. Brayden Baby has been designed to be an ideal teaching manikin for Paediatric BLS, EPALS and EPILS. Developed with help from leading European experts in the field of infant resuscitation, Brayden Baby CPR Manikin's unique features include intuitive lights that provide real-time feedback for ventilation training and lights providing real-time feedback for chest compressions (depth, rate, recoil and finger/thumb position). When all component parts of chest compressions are performed correctly, a CPR quality indication light illuminates on the forehead of the Brayden Baby manikin. Both sets of lights help guide the student to perform good quality infant CPR. **Aero Healthcare is the exclusive Australian importer of the entire [Brayden range](#), including the Brayden Baby, the [Brayden Pro](#), [Brayden Pulse](#), [Brayden Junior](#) and the [Brayden](#).**



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes