



## LAERDAL QCPR Manikin Baby Anne

\$667.39

In stock  
SKU LMQ-B

 Add to Wish List

 Add to Compare

LAERDAL QCPR Manikin Baby Anne



### Details

For years, instructors of infant CPR have not been able to know for sure what the participants learned. Until now. Little Baby QCPR is a realistic and affordable BLS manikin for pediatrics. By utilizing gamified learning and objective feedback, Little Baby QCPR improves training quality, learner engagement, and classroom efficiency.

### Easy to use. Quality assured

**Set up a class in no-time** Use the one-click connection to set up training. The QCPR app onboarding guides you through the setup of your classroom. **Engage your learners** Increase engagement and confidence in both learner and instructor with a combination of realistic feel and intuitive feedback. **Guideline compliant** With objective feedback on compressions and ventilations, Little Baby QCPR is future-proof. It teaches all the parameters of high-quality CPR as defined by the AHA.

### Install confidence in your learners

**Foreign body removal** Choking is one of the biggest reasons for cardiac arrest in infants. With Little Baby QCPR you allow learners to practice the not so gentle touch to see what it takes to remove airway obstructions. **Correct hand placement**



---

Spot incorrect hand placement with built-in sensors and QCPR app to instantly adjust learners performance. **Act when it matters** Practice correct compressions, ventilations, and foreign body removal so that learners act with confidence when it matters the most.



### More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes