



AERODIAGNOSTIC Digital Clinical Thermometer

\$7.60

In stock
SKU MT80D

 Add to Wish List

 Add to Compare


Digital Clinical Thermometer



Details

Size: 13.5cm **Pricing Unit:** unit/1 **Carton QTY:** 240 The AERODIAGNOSTIC range offers a selection of commonly used diagnostic supplies in the healthcare industry.

Thermometers

Thermometers are critical for measuring body temperature to diagnose and monitor fever or illness. Available in digital or glass thermometers. 



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes