



AERODIAGNOSTIC Diagnostic Pen Light



\$3.81

In stock
SKU PTD06

 Add to Wish List

 Add to Compare


Disposable Penlight



Details

Size: 10cm **Pricing Unit:** unit/1 **Carton QTY:** 360 The AERODIAGNOSTIC range offers a selection of commonly used diagnostic supplies in the healthcare industry.

Diagnostic Penlights

Diagnostic penlights are useful for providing a light source for testing pupil dilation and also as a secondary light source when examining first aid injuries. 



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes