




## QUIKCLOT Haemostatic EMS Roll 7.5cm x 1.2M

**\$70.73**

In stock  
**SKU** QC-EMSR

 Add to Wish List

 Add to Compare  
QuikClot EMS Roll



### Details

QuikClotEMS Rolled Gauze is 3 inches x 4 feet of soft, white, nonwoven, hydrophilic gauze impregnated with kaolin, an inorganic mineral that is both safe and effective in accelerating the bodys natural clotting cascade without any exothermic reactions or use of animal or human proteins. Each QuikClotEMS Rolled Gauze, is packaged in an easy-tear pouch.

- Effective, Simple, and Safe
- Easy to use and remove
- Flexible and pliable

The QuikClot range of hemostatic products have been chosen by EMS, Law Enforcement, Military and Healthcare professionals to help save lives everyday.



### More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes