



## AERORESCUE Basket Aviation Stretcher

**\$1,849.83**

In stock  
**SKU** STR-08

 Add to Wish List

 Add to Compare

Basket Aviation Stretcher



### Details

**Size:** 220 x 64.5 x 16.5cm **Pricing Unit:** unit/1 **Carton QTY:** 1 AeroRescue Stretchers are a light-weight high strength solution for many and varied patient recovery situations. Available in a variety of designs and materials for use in a wide range of environments and applications.

### STR-08

**Weight limit:** 205kg Sturdy and flexible with quick reliable fittings and special sling equipment, the basket aviation stretcher is ideal for all types of aid in special situations such as lifting and transport by helicopter. Materials are fireproof and highly resistant to wear and corrosion.



### More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes