



## Wooden Non-Sterile Tongue Depressors Box/100



**\$6.78**

In stock  
**SKU** TDNS001

 Add to Wish List

 Add to Compare

Wooden Tongue Depressors



### Details

**Size:** 15cm **Pricing Unit:** box/100 **Carton QTY:** 50 The AERODIAGNOSTIC range offers a selection of commonly used diagnostic supplies in the healthcare industry.

### Tongue Depressor

A wooden tongue depressor is used to depress the tongue when examining the throat or mouth.



## More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes