





AEROINSTRUMENTS Stainless Steel Universal Shears with Plastic Tip 15cm

\$2.87

In stock
SKU USP15-P

 [Add to Wish List](#)
 [Add to Compare](#)
Universal Scissors



Details

Size: 15cm **Pricing Unit:** unit/1 **Carton QTY:** 480 The AEROINSTRUMENT range offers a wide range of surgical instruments manufactured from high-quality stainless steel components to suit first aid requirements. Assorted sizes and specifications have been selected for use in common first aid procedures.



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes