



QV Gentle Conditioner - 500g

\$21.33

In stock
Only **499** left
SKU 11728

 [Add to Wish List](#)

 [Add to Compare](#)

QV Hair Gentle Conditioner is designed to care for those that experience scalp sensitivity. The unscented conditioner's pH-balanced formulation is dermatologically tested and gently infuses moisture to the hair.



Details

QV Hair Gentle Conditioner is designed to care for those that experience scalp sensitivity. The unscented conditioner's pH-balanced formulation is dermatologically tested and gently infuses moisture to the hair.

Directions:

After cleansing with QV Hair Gentle Shampoo, massage QV Hair Gentle Conditioner evenly through the hair, leave for 1 minute, then rinse well. For more intensive conditioning, leave your hair conditioner in for at least 5 minutes before rinsing.

Benefits:

- Free from fragrance.
- Free from soap and colour.



- pH balanced.
- Ideal for sensitive skin and dry scalps.
- Vegan-friendly - contains no animal-derived ingredients.

Caution:

Avoid eyes. If product enters eyes flush with water. Store below 30°C.

Ingredients:

Aqua (Water), Glycerin, Cetearyl Alcohol, Panthenol, Behentrimonium Methosulfate, 1,2-Hexanediol, Caprylyl Glycol, Titanium Dioxide, Citric Acid.