



ZOLL AED 3 BLS Semi-automatic Defibrillator

\$4,693.92

In stock
SKU
8532-001203-13

 [Add to Wish List](#)

 [Add to Compare](#)

Professional-grade AED for first responders, with a full-colour ECG display, the Real CPR Help dashboard and one of the fastest shock-after-compression times available.

Details

The ZOLL AED 3 BLS is built for professional and first-responder use. It guides the rescuer through high-quality CPR and is one of the fastest AEDs at delivering a shock once compressions stop — two things that matter most to a person's chance of survival.

Best for: ambulance services, fire and rescue, surf and ski patrol, event medical teams, and trained workplace responders.

Includes:

- Large colour touch-screen display showing the patient's ECG
- Enhanced Real CPR Help with a CPR dashboard: compression rate and depth, a CPR cycle countdown timer, shock count and total elapsed time
- RapidShock Analysis for one of the shortest rhythm-analysis pauses available
- CPR Uni-pads universal adult and paediatric electrodes (5-year shelf life)
- Smart Battery Pack (5-year shelf life)
- ZOLL AED 3 carry case plus the new soft carry case
- Wi-Fi enabled for PlusTrac cloud readiness monitoring and Case Review
- One free year of PlusTrac, operator's manual, quick-start guide and rescue poster
- 8-year standard warranty

As an authorised Zoll supplier we also stock the matching pads, batteries, cabinets and PlusTrac program management.