



## ZOLL AED 3 BLS Aviation Defibrillator

**\$4,693.92**

In stock  
**SKU**  
8562-001206-13

 [Add to Wish List](#)

 [Add to Compare](#)

The aviation specification of the ZOLL AED 3 BLS, built around an FAA-certified battery for commercial and charter aircraft fit-out.

### Details

The ZOLL AED 3 BLS Aviation model is the professional BLS unit prepared for the cabin environment. It pairs the same colour ECG display, Real CPR Help dashboard and RapidShock Analysis with an FAA TSO-C142b-certified battery, tested to withstand vibration, altitude, temperature and EMI for reliable fleet readiness.

**Best for:** airlines, charter and aero-medical operators fitting out aircraft to aviation safety requirements.

#### Includes the BLS feature set:

- Large colour touch-screen display showing the patient's ECG
- Enhanced Real CPR Help with the CPR dashboard (rate, depth, cycle timer, shock count, elapsed time)
- RapidShock Analysis for one of the shortest rhythm-analysis pauses available
- CPR Uni-padz universal adult and paediatric electrodes (5-year shelf life)
- FAA TSO-C142b-certified aviation battery pack in place of the standard pack
- Wi-Fi enabled for PlusTrac cloud readiness monitoring and Case Review
- One free year of PlusTrac, operator's manual, quick-start guide and rescue poster
- 8-year standard warranty

Talk to us about fleet quantities and PlusTrac program management. As an authorised Zoll supplier we stock the matching consumables.