



## ZOLL AED 3 Battery Pack

**\$277.37**

In stock

**SKU**

8000-000696

 Add to Wish List

 Add to Compare

The genuine smart battery for the ZOLL AED 3, with a 5-year standby life and self-test reporting through the device.



### Details

The genuine ZOLL AED 3 battery pack keeps the device rescue-ready and reports its own status through the AED 3 self-test and, where fitted, PlusTrac readiness monitoring. We recommend keeping a spare for any AED held in a high-readiness setting.

**Best for:** AED 3 owners replacing an expired or used battery, or holding a spare.

Includes one ZOLL AED 3 battery pack, 5-year shelf life. Compatible with the standard ZOLL AED 3 (not the Aviation model).