



## AeroPlast Blue Dressing Length 7.5cm x 1m Xray detectable

**\$30.67**

In stock  
Only **1000** left  
**SKU 9005L**

 Add to Wish List

 Add to Compare

### Details

- Hypo-allergenic
- Latex Free
- Waterproof
- Metal Detectable
- Sterile

Any form of food handling demands that the strictest hygiene regulations are enforced. Contamination of any scale entering the food chain is potentially disastrous and extremely costly. AeroPlast Detectable provides a cost effective wound care solution with multiple forms of detection safe guarding your business against avoidable litigation.



Ideal for Manufacturing and Processing Sector.

- Bright Blue
- Metal Detection
- High Level of Adhesion
- Comformable
- Hypoallergenic Adhesive (Latex Free)
- Waterproof
- Breathable
- Centralised Wound Pad
- Sterile

**Dimensions:** 7.5cm x 1m

### **Important**

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

### **Related Products**

Check items to add to the cart or select all