





Blood Pressure Monitor Economy with Adult Cuff 2 Tubes

\$91.26

In stock
Only **1000** left
SKU 400310

 [Add to Wish List](#)
 [Add to Compare](#)

Details

The Economy blood pressure monitor features:

- Modern design and a 300mmHG manometer with no pin stop mechanism
- Nylon fabric cuff with velcro
- Chrome plated inflation valve
- Sphygmomanometers supplied with zippered leatherette storage case

Note: This unit comes without a stethoscope included.

Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

Related Products



Check items to add to the cart or select all