



## BVM Silicone Resus Set Rescu-2 Child (2 Masks + 2 Airways) (LRI)

**\$241.02**

In stock  
Only **1000** left  
**SKU** 440011

 Add to Wish List  
 Add to Compare

### Details

Rescu-2 is an economical durable silicone manual resuscitator which has been designed to meet the need for high performance. It is constructed from polycarbonate and silicon material and is autoclavable up to 134°C.

Features:

- Child resuscitator 550ml with 40cm H2O pop-off valve
- 2500ml oxygen reservoir
- 2 child masks #5133 and #5122
- 2 guedel airways
- See through carry case with instruction manual

### Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

### Related Products

Check items to add to the cart or select all