



## Chill-Its Evaporative Cooling Bandana with Cooling Towel - Solid Blue

**\$15.86**

In stock  
Only **1000** left  
**SKU** 120056

 Add to Wish List

 Add to Compare

### Details

Evaporative Cooling Bandana with Cooling Towel: Advanced PVA material keeps you cool without the bulk.

Construction, Trades, Maintenance, Landscaping/Grounds, Assembly/Fabrication, Material Handling, Freight/Baggage, Warehousing/Distribution, Iron/Steel Fabrication

- Advanced PVA material activates easily and quickly when run under water for one minute
- Re-usable: just soak in water to re-activate
- Super-evaporative cooling process helps keep you cool
- Lightweight and low profile; perfect for under hard hats
- Tie closure
- Patented design

### Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

### Related Products

Check items to add to the cart or select all