



## Clear Plastic Biohazard Waste Bags 260 x 160mm with Red writing

**\$2.39**

In stock  
Only **1000** left  
**SKU** 99236

 Add to Wish List  
 Add to Compare

### Details

Biohazard Plastic Biohazard Bags consists of two compartments - a zip-lock sealed section for storing specimens, test tubes, swabs, etc and an unsealed/open section for message inserts (i.e. request for testing or analysis). Bags are clearly marked with biohazard symbol for easy identification and handling.

**Size:** 260 x 160mm

### Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury

### Related Products

Check items to add to the cart or select all