



## Emesis Vomit Bags (50) with Plastic Bowl

**\$50.66**

In stock  
Only **999** left  
**SKU** 99702

 Add to Wish List

 Add to Compare





## Details

- Sentry Medical Emesis Vomit Bag has a firm polyethylene twist and lock rim to provide easy control and help prevent spillage after use. They isolate splash back and have a graduated indicative measurement system to 1500ml, allowing measurement of contents in the bag. On completion of use, the bag is easily disposed of with no mess. Simply twist the bag and lock into the keyhole slot on the rim. Rims are colour coded red for infection control.

### Contents:

- Pack of 50
- Plastic bowl

### Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.