



Intensity 5000 Hybrid Tens Machine

\$191.94

Out of stock
SKU CST5000

 Add to Wish List

 Add to Compare



Details

The InTENSity 5000 Hybrid TENS Machine is both digital and analog. This device combines the ease of traditional knobs with the beauty of LED technology to give the user the feedback they need to get the right therapy at the right time. The InTENSity™ 5000 Hybrid has 5 modes, a timer and comes in a hard case.

The electronics of the InTENSity 5000 create electrical impulses whose intensity, duration, number per second and modulation may be altered with the controls/switches. These buttons are very easy to use and the LED display showing the exact mode, time and intensity.

Includes:

- 1 x InTENSity 5000 unit
- 2 x Electrodes Leads
- 4 x 40mm x 40mm adhesive electrodes
- 1 x 9V Alkaline Battery, type 6F22
- 1 x 240v AC Adapter
- 1 x Instruction Manual
- 1 x Carrying case



Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

Related Products

Check items to add to the cart or select all