



## Laryngoscope Set Liberty Macintosh Style Blades (4) and handle

**\$596.77**

In stock  
Only **997** left  
**SKU** 400173-  
BYDAD0201

 [Add to Wish List](#)

 [Add to Compare](#)

### Details

Stainless steel, with sure-grip handle, the Liberty fibre optic laryngoscopes are high quality and autoclavable. Super fine quality, with Interchangeable stainless steel blades, fibre optic laryngoscope set Macintosh type.

### Features:

- Stainless steel
- Autoclavable to 134B0C
- Cold fibre optic light
- Macintosh Blades in four (4) sizes: Infant; Child; Adult and Large Adult
- Fully compatible with all fibre optic hook-on handles that conform to ASTM (F965)
- Halogen bulb
- Sure grip handle



- Satin non-glare finish
- Takes 2 x C size batteries (not included)

### **Important**

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

### **Related Products**

Check items to add to the cart or select all