



Nitrile Finger Cots - Anti-static/Powder free - Blue - 100pcs - Large

£20.72

In stock
Only **1000** left
SKU 230036

 Add to Wish List

 Add to Compare



Details

- Nitrile - Latex Free and blue in colour for easy split detection
- Pre - rolled for convenience, comfort and protection
- Durable and long wearing for cost effectiveness
- They are ideal for restaurants, cafés, butchers, food manufacturing and all food preparation.

Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

Related Products

Check items to add to the cart or select all