



## Opsite Flexigrid Transparent Film Dressing - Singles

**£1.08**

In stock  
Only **985** left  
**SKU** TM-SN4628

 Add to Wish List  
 Add to Compare

### Details

Dressing Moisture Vapour Permeable Adhesive Film Dressing.

OpSite Flexigrid is a transparent, adhesive film, with a unique wound measurement grid. The film is moisture vapour permeable, conformable and extensible. It is clinically proven to provide moist wound healing for superficial wounds; it is also ideal for use as a secondary dressing and to help prevent pressure sores.

### Features:

- Faster, safer, less painful wound healing
- Superior film with polyacrylate adhesive
- Waterproof film
- Simple, effective application technique
- Adaptable to awkward areas Unique wound measurement grid

### Benefits:

#### **Faster, safer, less painful wound healing**



- OpSite Flexigrid has set the standard in moist wound healing by retaining wound contact with the natural wound exudate which contains vital nutrients for growing cells, and white blood cells to prevent infection. OpSite Flexigrid provides an environment for safer and less painful wound healing. The high moisture vapour permeability allows excess exudate to evaporate, preventing skin maceration

### **Superior film with acrylic adhesive**

- The OpSite Flexigrid film is highly extensible and conformable, which combined with its high moisture vapour permeability, means the dressing can be left in place for up to 14 days. This minimises interference at the wound site and reduces the amount of nursing time required

### **Waterproof**

- OpSite Flexigrid is impermeable to water and body fluids, protecting the skin underneath and allowing the patient to bathe without changing dressings. It also aids in the prevention of bacterial contamination

### **Simple, effective application technique**

- OpSite Flexigrid provides a simple and effective application system. The flexible carrier enables single handed application allowing the dressing to conform to the contours of the body. Its transparency enables accurate positioning over the wound site

### **Adaptable to Awkward Areas**

- OpSite Flexigrid has been designed to be adaptable for use on all areas of the body. It can be easily cut into strips or shapes for application to awkward areas such as heels, elbows and knees

### **Unique Wound Management Grid**

- OpSite Flexigrid incorporates a unique wound measurement device which allows simple but accurate assessments of the healing process to be made. After tracing the edges of the wound, the OpSite Flexigrid carrier can be removed and stored in the patient's notes

### **Options: (Singles)**

6 x 7cm

10 x 12cm

### **Important**

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.



---

**Related Products**

Check items to add to the cart or select all